

PRSD-SERIES

Precision Rotation Stages. 360 degrees coarse adjustment and fine adjustment within 10 degrees with centre aperture allowing light beams to go through. The best applications for this stage are polarization and image adjustment.

ITEM	PLATFORM SIZE	ACTUATORS	CENTRE APERTURE	AXIS OF TRAVEL	ECCENTRICITY	COARSE ADJUSTMENT	FINE ADJUSTMENT	VERNIER SCALE	WEIGHT
PRSD-80A	Φ 80mm	100TPI	Φ=40.0mm	Θ	25um	360°	5°	0.1°	
PRSD-100A	Φ 100mm	100TPI	Φ=50.8mm	Θ	25um	360°	5°	0.1°	

PRSD-80A



PRSD-100A

